

# Development of Next-Gen Smart Mat to Measure Oxygenation and Ulcer Risk

Jeffrey Schoess<sup>1</sup>, David G. Armstrong<sup>2</sup>, Karen D'Huyvetter<sup>2</sup>  
Eden Medical, Inc.<sup>1</sup>, Howard Lake, MN., Keck School of Medicine of USC<sup>2</sup>

## The Problem



### Common Pathway of Diabetic Foot Ulcer Occurrence and Recurrence



#### Key Facts-

- 9.1 to 26.1 million diabetic foot ulcers develop per year worldwide<sup>1</sup>
- 785 million outpatient visits (between 2007 and 2013) due to DFU and infection<sup>2</sup>
- \$176 billion per year in US for diabetic direct costs<sup>3</sup>
- 40% DFU recurrence within 1 year (60% within 3 years)<sup>4</sup>
- Strongest predictor of DFU is previous ulcer<sup>5</sup>
- Key Factors- early signs of skin damage (callus, blistering, or hemorrhage), repetitive shear stress, loss of sensation
- Early recognition is a key to preventive care- oxyhemoglobin (H<sub>b</sub>O<sub>2</sub>) and tissue oxygen saturation (S<sub>t</sub>O<sub>2</sub>) to characterize microcirculation dysfunction and oxygenation changes, the early biomarkers of ulcer formation.


## Smart Mat Approach (730/850 nm)



$$A = \epsilon c l$$


- A = absorbance (no units).
- $\epsilon$  = Molar extinction coefficient (M<sup>-1</sup> cm<sup>-1</sup>).
- c = Concentration (M).
- l = pathlength (cm).

#### Beer-Lambert Law



Light Pipe Pic

## Smart Mat




Smart Mat Design Layout

## Smart Mat Storyboard


### Concept ID

#### Sensor Array




Occlusion Test  
(58 year old Type 2 diabetic)

## Bench Test Results



HbO2 Plot



Hb Plot

## Lessons Learned

- Dual Wavelength Near Infrared Method- Is sensitive enough to quantify and track the early biomarkers of ulcer formation (i.e. oxyhemoglobin, deoxyhemoglobin)
- Wide Dynamic Range- our dual wavelength NIRS technology has a dynamic range of 96 db (i.e. 16 bit resolution capability)
- System Responsivity- the input/output gain of the optical system is adequate to detect the hemoglobin markers for a diabetic subject

## Conclusions

- Foot Ulcer Occurrence/Recurrence- is common and costly
- There is room for improvement in the effective management of diabetic foot ulcers, to reduce the risk of re-ulceration
- The Smart Mat is an innovative and new approach for reducing ulceration risk in the home setting to increase ulcer-free days

## References

- Diabetes atlas. 7th ed. Brussels: International Diabetes Federation, 2015. www.diabetesatlas.org)
- Skrepnek GH, Mills JL Sr, Lavery LA, Armstrong DG. Health Care Service and Outcomes Among an Estimated 6.7 Million
- Ambulatory Care Diabetic Foot Cases in the U.S. Diabetes Care 2017 May 11
- Armstrong DG, Boulton AJ, Bus SA, Diabetic Foot Ulcers and their Recurrence, N Engl J Med. 2017 Jun 15;376(24):2367-2375